



# Titan power The best just got better

---

Impedance, OAE, ABRIS  
and the revolutionary  
Wideband Tympanometry  
in one powerful device



**Interacoustics**

# Let your equipment and needs interact

The first and only one of its kind 4 superior modules in 1 device

---

**The Titan from Interacoustics is a modular platform that lets you combine the tools you need with handheld freedom and full flexibility - from one trusted supplier.**

---

## One size fits all

- Screening
- Diagnostic
- Advanced Clinical

Customize your Titan for screening, diagnostic and advanced clinical testing. Match the modules with your current needs - and upgrade when your challenges change. The Titan delivers Impedance, ABRIS, OAE and the revolutionary Wideband Tympanometry.

---

## Built for the future

Every Titan software module is developed based on proven technologies and feedback from users and audiology experts. Features are developed in sync with your changing needs and contain the latest technologies making the Titan a future-safe investment.



---

## A preview of the benefits

The Titan is truly a mighty platform with ease of use, professionalism and convenience in one device. Look forward to:

- Superior, professional aesthetics
- Lightweight
- High resolution color display
- Patient-centered testing through multiple probe and transducer configurations
- Handheld and PC mode for full flexibility
- Easy to use with even a minimum of technical understanding
- Supreme usability through storage of calibration values in probe and ID-transducers



---

**Impedance, OAE,  
ABRIS & Wideband  
Tympanometry in  
one powerful device**

**Status lights**  
Multi colored  
LEDs display test  
status for  
ease of use.



# Set your winning combination on one mighty platform

Customize your Titan through dedicated software modules to create a winning combination for screening, diagnostic and advanced clinical testing.

## **IMP440**

Screening, diagnostic and clinical impedance

## **ABRIS440**

Automated ABR for infant screening

## **DPOAE440**

Screening and clinical OAE

## **TEOAE440**

Screening and clinical OAE

A person wearing a white lab coat is shown from the waist up. A handheld device, the Titan, is tucked into the right pocket of the lab coat. The device is black and white with a blue top. The background is a dark, textured surface on the left and a light grey background on the right. A white circular callout box is positioned in the upper right quadrant, containing text.

**Easy handling**

From standard screening to advanced diagnostics, the Titan is created with superior usability in mind. From pre-set protocols to customized settings, you will benefit from intuitive interfaces and a professional design.



Pre loaded or  
customized?  
Use existing or  
create your own  
protocols and  
print templates

# IMP440 for every thinkable task

## Your multi-tasking business partner

The IMP440 is available in a screening, diagnostic or clinical version. It features standard tympanometry, ipsilateral and contralateral acoustic reflexes, reflex decay and reflex latency, as well as three Eustachian tube function tests.

## The fastest, most reproducible tympanometry available

IMP440 delivers tympanometry through advanced technologies. These include: "Endless airflow" pressure control technology, high resolution recordings and an intelligent pump system with an adaptive speed control around the tympanic peak. The IMP440 module will help streamline and improve your daily workflow - in record speed.

## IMP440 comes with...

**Three Eustachian tube function tests**  
Intact, perforated and patulous

**Ipsi and contra acoustic reflexes**  
Quick screenings at a single intensity or fast threshold searches

**Reflex decay**  
Ipsilateral as well as contralateral stimulation


**Reflex latency**  
Detailed investigation at the onset of the reflex

**Optional high frequency probe tones**  
(678, 800 and 1000Hz)

**Shoulder unit and probe indicator LED**  
For ease of use and visual feedback of status

## Print wizard

The embedded print wizard in the Titan suite can merge data from all four modules and other Interacoustic Suites for full patient reporting.



---

# ABRIS440 Automated ABR in seconds

---

## Automated ABR

The ABRIS440 module for Titan is fast and easy automated ABR testing for infants and adults. The screening module uses the ground breaking CE-Chirp® stimulus and weighted averaging to reduce test time by up to 50% compared to traditional click stimuli. With an easy to understand display (PC or handheld device), HiTrack compatibility and user customizable protocols, the ABRIS440 combines fast, quality testing with unrivaled usability.

---

## Simply press "START"

Following patient preparation, simply press start and await a pass/refer result. Minimal test time is required and the ABRIS440 default CE-Chirp® (35dBnHL) protocol provides a reliable automated ABR test with a sensitivity of 99.9% and a specificity of more than 96%. Choose from four transducers that are automatically detected when connected to the preamplifier.

### CE-Chirp®

Stimulus compensates for cochlear travel time. The result is a response amplitude up to twice the size compared to a standard click.



# DPOAE440 TEOAE440

## Easy screening and diagnostics

**OAE modules with specialized advantages: diagnostic OAE assessment and newborn hearing screening.**

### With you through any challenge

The screening versions offer speed, precision and user-friendliness. The Titan ensures exact intensity stimulation throughout the test and monitors the ongoing noise resulting in the most accurate results achievable. Combined with Bayesian weighted averaging that counteracts excess patient noise, the software is extremely efficient especially for infant screening.

The clinical versions provide unsurpassed flexibility and a comprehensive frequency range for a full diagnostic evaluation. Customizable protocols offer an unlimited number of tests utilizing a multitude of parameters including stimuli, intensity, stop criteria and much more.

Both modules are capable of performing pressurized OAE measurements, which is especially useful when testing infants & younger children with middle ear conditions. The technique reduces the effects of negative middle ear pressure and more children can potentially be diagnosed earlier.

Below:  
DPOAE440  
DP-Gram  
with a PASS  
in bar view



### Screening highlights

#### DPOAE440 & TEOAE440

- DP range 0.5-6kHz
- TE range 1-4.5kHz (60-84 dB SPL)
- Customizable PASS/REFER criteria
- Advanced AGC
- Bayesian weighted averaging
- Noah & Hitrack compatibility
- 226Hz tympanometry
- Operate from Titan or PC (optional)
- Online quick guides for fast learning
- Pressurized OAE (optional)

### Clinical highlights

#### DPOAE440 & TEOAE440

- DP range 500Hz-10kHz
- TE range 500Hz-5.5kHz (30-90 dB SPL)
- Fully customizable protocols
- Basic & Advanced view
- Manually controlled test time
- Unlimited test points
- Advanced AGC
- Bayesian weighted averaging
- Noah & Hitrack compatibility
- DP Input/Output test
- 226Hz tympanometry
- Operate from Titan or PC
- Pressurized OAE



### Flexible freedom

The flexible design ensures that any screening program or advanced audiology department can adapt protocols to meet their exact test requirements. As a hybrid, the Titan can be operated either handheld or from your PC.


# Introducing a revolution Wideband Tympanometry

---

Wideband Tympanometry  
- WBT - expands our understanding of the middle ear.

WBT adds a new dimension in differential diagnostics.

---



---

### One sweep: Multiple tympanograms

Differentiate pathologies with the aid of a powerful 3D image of tympanometric data across a full range of frequency and pressure. The power of a click alternative to measure at peak pressure.

---

### The power of a click

The Click stimulus unlocks the power of Multi-Frequency Tympanometry in the same time as testing a single frequency.

The 3D result reveals information of middle ear performance at any frequency across a full range of pressure change.

---

### Combine with OAE and automated ABR

WBT provides missing data about middle ear conductive issues that will enhance a diagnosis based either on OAE or ABR alone.

---

### 3D graphs

Absorbance is revealed in a multi-dimensional graph that enables the clinician to have an immediate understanding of middle ear performance across frequency and pressure. Diagnostic interpretation is made easy with a clear map of middle ear function.


---

### No extra skills required

WBT is easy and yet comprehensive. Differential diagnostics are made much more intuitive with the aid of overlaid normative data. Compare to illustrated examples of elusive pathology.

**A clear benefit**  
WBT streamlines middle ear diagnostics and provides accurate and specific measures across a comprehensive frequency range.





## Interacoustics

Interacoustics is a world leading provider of diagnostic solutions in the field of hearing and balance assessment. Since 1967 we have designed and manufactured our innovative diagnostic solutions for the world of audiology with a constant focus on providing our customers with quality, dependable products.

---

## With you at all times

We operate in over 100 countries worldwide through a carefully selected network of distributors and service centres. Purchasing an Interacoustics product guarantees not only a quality product, but also direct access to our highly professional training and support service.

---

## Product specifications

All technical and hardware specifications concerning software modules and the Titan platform can be downloaded from our website.

Read more at  
[interacoustics.com](http://interacoustics.com)


---

## Interacoustics A/S

Audiometer Allé 1  
5500 Middelfart  
Denmark

T +45 6371 3555  
F +45 6371 3522

[info@interacoustics.com](mailto:info@interacoustics.com)  
[interacoustics.com](http://interacoustics.com)



Go online to  
explore our  
full product  
range



# Interacoustics